



Cat. No: 7.2.3 - 32 / 06 / June 2020

PROTECTIVE environment EASY to move

The TINC 101 is a lightweight, easily transportable incubator. Only two care providers are needed to handle it. The flexible design of the TINC 101 makes it suitable for all ambulances.

This incubator may be used in conjunction with ventilators and other equipment. It is provided with iris ports that guide all tubing including ventilation tubing. The wide ports of the acrylic infant hood permit convenient access to the infant.

The moulded, single-piece infant bed and the interior of the incubator are crevice free and very easy to clean. The bed may be slid readily to provide access to the head of the infant for swift intubation or resuscitation.

Multi-position restraining straps hold the infant securely during transport.

The TINC 101 features an optional electrically operated suction module and a user friendly control panel that offers strain-free

Upright

viewing. Other features of this incubator include a separate infant temperature probe and a convenient IV pole.


The incubator accepts a 5 litre oxygen cylinder and an auxiliary 5 litre cylinder. A quick-fitting pin mechanism connects to the cylinder.

A gooseneck lamp is provided to facilitate infant observation even under conditions of darkness. The lamp can provide illumination on any desired area.

The TINC 101 operates on its own battery, which can support it for 3 hours, or on power from a vehicle or AC mains. Its stand is optionally fixed or collapsible. With the collapsible option, the stand self-folds at its station in the transport vehicle.

Noise and vibration inside the infant crib are minimal.

The TINC 101 is backed by the many years of experience PHOENIX Has in designing and supplying sophisticated neonatal care equipment worldwide.



endeavours to meet your every need. If your practice has any special requirement, call the marketing department at PHOENIX. They will be happy to provide you a solution.

Phoenix Medical Systems Private Limited